

HAFFMANS IAM INPACK 2000 AIR METER

GENERAL PRODUCT INFORMATION

In the beer and beverage industries, the content of dissolved oxygen (O_2) is a determining factor in the quality, taste and flavor stability of beer and carbonated beverages, as it causes a rapid decline of the flavor stability and shelf life. By measuring the air volume in the container, the O_2 content can be determined. Breweries as well as soft drink manufacturers continuously control and measure the concentration of air after filling.

The Inpack 2000 Air Meter, type IAM, is used to determine the air content in the headspace volume as well as the total air content in the beverage, in combination with the Inpack 2000 CO_2 Calculator or the Inpack 2000 CO_2 Meter. The principle is based on CO_2 absorption by caustic.

The pressure of the container is released and the gas fed to the IAM. The CO_2 is completely absorbed by the caustic and the remaining gas is collected in the measuring burette. Its volume (ml) is shown on the measuring burette.

The device is available in two executions:

- Type IAM-20 has a measuring range of 0-20 ml, with higher accuracy in the low range and is typically used for carbonated beverages with low air content such as beer
- Type IAM-36 measures the range of 0-36 ml and is typically used for carbonated beverages with a higher air content, such as soft drinks



CUSTOMER BENEFITS

- **Perfect quality management**
- **Cost saving**
 - maintenance free

APPLICATIONS

- **Laboratory, in combination with Pentair Haffmans' Inpack 2000 CO_2 Calculator and Inpack 2000 CO_2 Meter**

HAFFMANS IAM

INPACK 2000 AIR METER

TECHNICAL DATA

IAM-20

Measuring range

0.0 – 20.0 ml

Subdivision

Range 1	0.0 – 0.50 ml: 0.05 ml
Range 2	0.6 – 6.50 ml: 0.10 ml
Range 3	7.0 – 20.00 ml: 0.50 ml

Accuracy

± 5 % of measured value with a min. of 0.025 ml

Caustic content

Approx. 160 ml

Weight

Ca. 1.4 kg

IAM-36

Measuring range

0.0 – 36.0 ml

Subdivision

Range 1	0.0 – 6.0 ml: 0.10 ml
Range 2	7.0 – 36.0 ml: 0.50 ml

Accuracy

± 5 % of measured value with a min. of 0.05 ml

Caustic content

Approx. 175 ml

Weight

Ca. 1.3 kg

Dimensions

250 x 170 x 595 (L x W x H mm)

Caustic concentration

30 % g/l NaOH with approx. 1 % methanol

30 % g/l KOH with approx. 1 % methanol

SCOPE OF SUPPLY

- Inpack 2000 Air Meter*
 - Type IAM-20 or
 - Type IAM-36
- Measuring burette and leveling vessel with clamp
- Set of tubing
- Instruction manual

* When ordering, please state which instrument you require.

OPTIONS

- Certificate of measurement

**HAFFMANS B.V.**

P.O. BOX 3150, 5902 RD VENLO, NETHERLANDS WWW.FOODANDBEVERAGE.PENTAIR.COM

All Pentair trademarks and logos are owned by Pentair. All other brand or product names are trademarks or registered marks of their respective owners. Because we are continuously improving our products and services, Pentair reserves the right to change specifications without prior notice.

Pentair is an equal opportunity employer.

haffmans-inpack-2000-air-meter-2106-en ©2021 Pentair, All Rights Reserved.

